

















Accessories for 2.1027.0100

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1027.0100

Qty.	Order no.	Description	
1 pcs	1.1027.0010	Eco KF Titrator	 <p>The new Eco KF Titrator with integrated magnetic stirrer und touch-sensitive display is ideal for routine water content determination. Predefined methods enable a smooth and easy instrument start-up. In addition to GLP-compliant printouts on paper or as PDFs, the Eco KF Titrator also offers the option of connecting balances or of sending the determination data to a PC via PC / LIMS reports.</p>
1 pcs	1.1029.0010	Solvent Pump	 <p>Solvent Pump is used with an Eco KF Titrator for reliably adding solvent and aspirating off titration vessel content without contact with the reagents.</p>
1 pcs	6.01600.200	Siphon Breaker	 <p>Prevents continued running of liquid out of the reagent bottle into the titration cell.</p>
1 pcs	6.01804.210	PVC tubing 4/6 0.35 m	 <p>For tubing connections to the OMNIS Solvent Module or the Solvent Pump.</p>

1 pcs	6.01805.130	3/250 mm PTFE tubing with screw nipple	<p>Aspiration tubing with screw nipple to connect the reagent bottle to the Siphon Beaker.</p> 
1 pcs	6.02110.000	Eco KF Titrator / Solvent Pump connecting cable	<p>Cable for connecting a Solvent Pump with an Eco KF Titrator, used for energy supply</p> 
1 pcs	6.02712.000	Protective cover for Eco Dosimat / Eco Titrator	<p>Protective cover for Eco Dosimat / Eco Titrator display.</p> 
1 pcs	6.02801.000	Cotton for adsorber tube	<p>Wad of cotton (4 pieces), suitable for OMNIS 601406010 adsorber tube.</p> 
1 pcs	6.03002.210	10 mL cylinder unit for Eco Dosimat / Eco Titrator	<p>10 mL cylinder unit for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.</p> 
1 pcs	6.05610.050	KFT equipment for Eco KF Titrator	<p>Complete accessory set for volumetric Karl Fischer titration</p> 
1 pcs	6.1446.040	Threaded stopper / M6	

1 pcs	6.1446.080	Threaded stopper / M8	
1 pcs	6.1446.090	Stopper / B-14/15 / M8	
1 pcs	6.1543.110	Dosing tip / M8 thread Transfer tip with M8 thread	
1 pcs	6.1543.120	Aspiration tip / M8 thread Aspiration tip if 6.1805.200 siphoning tubing is used	
1 pcs	6.1546.040	Piston tool for Eco Dosimat/Eco Titrator Tool for removing and mounting the piston.	

2 pcs **6.1602.105** **Bottle attachment / GL 45 / green**

For direct dosing from reagent bottles with GL 45 thread.



1 pcs **6.1608.023** **Amber glass bottle / 1000 mL / GL 45**

For Exchange Units. Bottle for auxiliary solutions.



1 pcs **6.1608.030** **Round glass bottle / 1000 mL / GL 45**



1 pcs **6.1619.010** **Adsorber complete, for Exchange Unit**



1 pcs **6.1623.000** **Overflow protection for Ti Stand**

For 803 Ti Stand



2 pcs **6.1805.090** **FEP tubing / M6 / 31 cm**
With light and kink protection



1 pcs **6.1805.100** **FEP tubing / M6 / 40 cm**
With light and kink protection



2 pcs **6.1805.200** **PTFE tubing / M8 / 0.5 m**



1 pcs **6.1808.050** **Tubing adapter olive / M8 outer**
1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 pcs **6.1819.020** **Tube**
For M6 threaded holes. For 6.1602.1XX bottle attachments



1 pcs **6.1819.050** **Tube**
For M8 threaded holes of 6.1602.10X Bottle attachments with 6.1608.030 Aspiration bottle. For 703 and 803 Ti Stand, 681 Pump Unit and 679 Rancimat



1 pcs **6.2013.010** **Clamping ring**

For support rods with a diameter of 10 mm.



1 pcs **6.2016.070** **Support rod / 400 mm**



1 pcs **6.2023.020** **Clip for SGJ 14/15**
Clip for SGJ 14/16



1 pcs 6.2164.010 Power supply unit 100 - 240 V/24 V DC

Power supply unit for 915 KF Ti-Touch and 916 Ti-Touch



1 pcs 6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.










1 pcs 6.6081.020 USB flash drive for Eco KF Titrator

Storage medium for Eco KF Titrator



Optional accessories

Order no.	Description	
6.01807.000	Absorption cartridge	
	Absorption cartridge for OMNIS Solvent Module or Solvent Pump, for filling with 6.2811.000 molecular sieve.	
6.02709.040	Stopper	
	Stopper for titration vessel lids 6014050X0.	
6.02711.000	Spoon for paste	
	Suitable for use in Karl Fischer water content determinations with pasty substances.	
6.03002.150	5 mL cylinder unit for Eco Dosimat / Eco Titrator	
	5 mL cylinder unit for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.	
6.03002.220	20 mL cylinder unit for Eco Dosimat / Eco Titrator	
	20 mL cylinder unit for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.	
6.03002.250	50 mL cylinder unit for Eco Dosimat / Eco Titrator	
	50 mL cylinder unit for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.	
6.05610.060	KFT equipment for Eco KF Titrator with Polytron	
	Complete accessory set for volumetric Karl Fischer titration with Polytron	